

Automated STR Data Analysis: Validation Studies

Mark W. Perlin (Cybergenetics, Pittsburgh, PA)
David Coffman (Florida Department of Law Enforcement, Tallahassee, FL)
Cecelia A. Crouse & Felipe Konotop (Palm Beach County Sheriff's Office, FL)
Jeffrey D. Ban (Division of Forensic Science, Richmond, VA)

Automated Analysis
Databasing Validation
Casework Studies



Cybergenetics

NIJ grants 2000-IJ-CX-K005 & 2001-IJ-CX-K003

Reviewing STR Data

Human data review bottleneck

Computer Automation
Quality Assurance
Database Integrity
Casework & Mixtures

Key goals:

- no error
- high throughput
- small staff

TrueAllele™ Technology

Eliminates STR human review bottleneck

**Gel-based, or
Sequencer, or
Capillary**

**Raw STR
Data INPUT**

Fully Automated
(on Mac/PC/Unix)

Color Separation
Image Processing
Lane Tracking
Signal Analysis


Ladder Building
Peak Quantification
Allele Designation
Quality Checking
CODIS Reporting

**Quality
Assured
Profiles**

**Database
OUTPUT**

Protected by US patents 5,541,067 & 5,580,728 & 5,876,933 & 6,054,268

Multi-Platform Engine



Validation Methods

1. Obtain original data
2. Process data in TrueAllele ES
(auto-setup, process run, Q/A,
call alleles, apply rules, check)
computer: accept/reject/edit
3. Review all data
one person, many computers
human: accept/reject/edit
4. Generate results & stats

Rule Settings

		Gel		Capillary	
		On/Off	Value	On/Off	Value
Extract	Dispersion	On	60%	On	80%
	Extra Allele	On	25%	On	15%
	High Signal	On	9000	On	10000
	Low Signal	On	150	On	300
	Low Homozygote	On	250	On	500
	Third Peak	On	23%	On	15%
Amplify	Check Control	On	----	On	----
	Conflict	On	----	On	----
	High n Peak	On	30%	On	30%
	Negative	On	----	On	----
	Noise	On	----	On	----
	Relative Area	On	50%	On	50%
Separate	Relative Height	On	50%	On	50%
	Stutter	On	20%	On	15%
	Lane to Lane	Off	----	Off	----
	Off Ladder	On	0.4	On	0.4
Other	Peak Morphology	On	50%	On	20%
	Uncorrelated	On	0.4	On	0.4
	Amelo	On	----	On	----
	Crossover	On	----	On	----
Other	Dye to Dye	Off	----	Off	----
	Overlap	On	----	On	----
	Rare	On	----	On	----
	Unexpected	On	25%	On	15%

Hitachi + PowerPlex

Hitachi FM/Bio2 & Promega PowerPlex 1.2
~8,000 PBSO genotypes reviewed
TrueAllele performed all gel & allele processing

Computer: ~75%* data, no review needed
Human: All these designations correct

**TrueAllele expert system can eliminate
most human review of gel STR data**

Hitachi Results

Human Review

		Accept	Edit	Reject
Computer Process	Accept	72.5%*	0.0%	0.0%
	Edit	23.0%	4.2%	0.3%
	Reject	0.0%	0.0%	0.0%

non-automation data

ABI/310 + Pro/CoFiler

ABI 310 & ProfilerPlus/Cofiler
~24,000 FDLE genotypes reviewed
TrueAllele performed all CE & allele processing

Computer: ~85% data, no review needed
Human: All proper designations correct

**TrueAllele expert system can eliminate
most human review of CE STR data**

310 Results

Human Review

		Accept	Edit	Reject
Computer Process	Accept	86.4%	0.0%	0.0%
	Edit	10.9%	2.1%	0.3%
	Reject	0.3%	0.0%	0.0%

94 genos/min

ABI/3700 + Pro/CoFiler

ABI 3700 & ProfilerPlus/CoFiler
~17,000 FDLE genotypes reviewed
TrueAllele performed all CE & allele processing

Computer: ~85% data, no review needed
Human: All TA/ES designations correct

TrueAllele expert system can eliminate
most human review of CE-array STR data

3700 Results


Human Review

		Accept	Edit	Reject
Computer Process	Accept	84.9%	0.0%	0.0%
	Edit	10.8%	3.3%	0.8%
	Reject	0.2%	0.0%	0.0%


76 genos/min

The UK FSS Experience

Generate STR Data



TrueAllele expert system scores all STR data and assesses data quality



Person reviews a fraction of the data

UK National DNA Database

FSS ABI/377 Validation

Resources


- Data: 22,000 genotypes (SGMplus)
- People: 6 reviewers + 6 managers
- Time: 8 weeks work + 4 weeks report

Components

- Peak height correlation (GS vs TA)
- Establish baseline height (error-free)
- Designation accuracy (human vs TA)
- Network/computer environment
- QMS documentation

Results

- Greater yield with TA
- No errors on quality data



Casework Studies

Nonmixture
 Mixture
 Rape Kits
 Disasters
 LCN, SNPs

M.W. Perlin and B. Szabady, "Linear mixture analysis: a mathematical approach to resolving mixed DNA samples," *Journal of Forensic Sciences*, November, 2001.

Statistical Information

Prob{STR profile & peak quants}

Degenerate SGM+ alleles (6x6x6x10x...x6)

☑ 100,000,000 feasible profiles

Compute unknown minor contrib profiles

☑ @30% 1 feasible profile

☑ @10% 100 feasible profiles

LMA increases identification power
a million-fold; **CODIS match**

Validation Data Sets


Collaborators: Florida, Virginia,
New York, FBI, UK, Private Labs


Data: synthetic lab mixtures,
casework, rape kits, disasters


Input: TrueAllele quantitated peaks

Studies: comparison, concordance,
automated lab processes

Rape Kit: Unknown Suspect







Conclusions

- TrueAllele *databasing validation*
- Reduce time, error, staff & cost

- Ongoing *casework validation*
- Automate: data review & lab work
- Serve: police, courts, society
- Objective, comprehensive
